INSTRUCTIONS FOR USE

PROSTHETIC COMPONENTS

Medical Precision Implants S.A.
C/ Solidaridad, 2
P.I. Los Olivos
28906 Getafe - Madrid
España

(+ 34) 91 684 60 63
www.mpimplants.com

MPI Prosthetic components

Description:
Elements used for the fabrication of dental structures located on the implant.

MPI prosthetic components types:

MPI Privilege® internal connection (Codified as CM).
MPI Privilege® external connection (Codified as HE).
MPI Excellence® internal connection (Codified as EXC).

Cover screw or Healing screw

Medical device class IIb.
It is screwed directly to the implant to cap it.
Use: Temporary.
Types: CM and HE.
Material: Titanium.
Maximum torque recommended: 10 Ncm.

Healing Abutment

Medical device class IIb.
It is used in 2 phases surgeries to mould the gingiva. It is screwed directly to the implant.
Use: Temporary.
Types: Straight and aesthetic. CM, HE and EXC.
The aesthetic healing abutment serves to open more the gingiva and form it. It must be used with all the aesthetic prosthetic components so that the gum is not damaged.
Material: Titanium.
Maximum torque recommended: 10 Ncm.

Provisional

Medical device class I.
It is screwed directly to the implant to make a provisional prosthesis.
Use: Temporary.
Types: Hexagonal (Non-rotating) and cylindrical (rotating). CM and HE.
Material: material biocompatible PEEK. (not recommended for casting).
Maximum torque recommended: 10 Ncm.

Abutment screw / Gold screw

Medical device class IIb.
It is screwed directly to the implant to fix the prosthesis.
Use: Definitive.
Material (Abutment screw): Titanium
The final screw can be delivered sterilized (blister) to the dental clinics and not sterilized (bag) to the dental laboratories.
Material (Gold screw): Titanium with a gold plate.
The gold screw is always sterile (in blister).
Maximum torque recommended: 30 Ncm.

Implant Coping

Medical device class I.
It is screwed directly to the implant when it is integrated, for the coping taking.
Use: Temporary.
Types: HE (open tray), HE (closed tray), CM (open tray), CM (closed tray).
Open tray: there is a short screw to use it as a closed tray.
Closed tray: The post has 3 flat sides to facilitate its position in the silicone because it is a drag element.
Material: stainless steel.
Maximum torque recommended: 10 Ncm.

Replica

It is not considered as a medical device. It is used in the laboratory to simulate the implant in the model in order to make the prosthesis.
It is screwed to the implant coping.
Use: Temporary.
Types: CM and HE.
Material: stainless steel.

Conical abutment

Medical device class IIb.
It is used in the laboratory to make a cemented prosthesis.
It is screwed to the replica to create the prosthesis and then it is screwed to the implant.
Use: Definitive.
Types: Straight and angled (17° and 25°). CM and HE.
Material: Titanium.
Maximum torque recommended: 30 Ncm.

Cr-Co base Castable abutment

Medical device class IIb.
It is used in the laboratory to make a screwed prosthesis.
It is screwed to the replica.

Use: Definitive.

Types: Hexagonal (Non-rotating) and cylindrical (rotating). CM, HE and EXC.
Material: Cobalt-Chromium base and PEEK.
Maximum torque recommended: 30 Ncm.

Castable

It is not considered as a medical device. It is used in the laboratory to make a screwed prosthesis.
It is screwed to the replica.
Use: Definitive.
Hexagonal (Non-rotating) and cylindrical (rotating). CM and HE.
Material: Delrin.

Laboratory screw

It is not considered as a medical device. It does the function of the abutment screw in the laboratory and is delivered non-sterile.
It is screwed to the replica.
Use: Temporary.
Types: CM and HE.
Material: Titanium.
Maximum torque recommended: 30 Ncm.

Transethelial

Medical device class IIb.
It makes the function of a straight abutment to transform the internal connection into external connection or to elevate an infra-osseous implant connection.
It is screwed directly to the implant.
Use: Definitive.
Types: CM and HE.
Material: Titanium.
It has its own prosthesis. It is an intermediate abutment between the implant and the prosthesis.
For the uniblock straight abutment, the long screw M1,4 - TR009 is used.
For the angled abutment the short screw M1,4 - TR008 is used.
Our transepthelial system only allows structure restoration, because the castable of the system is only rotating.

Indications and intended use:

Prosthetic components are used to create a suprastructure on a dental implant to act as natural teeth.
Read the instructions carefully.
These devices are intended to be used in a sterile environment.
Personnel authorized to use these devices must have adequate knowledge, experience, education and training of these devices, accessories and surgical instruments.
Protective measures and necessary clothing must consider.
Check the expiration date written on the packaging and do not use if the sterile packaging is damaged or opened (in case of sterile devices).
Different prosthetic components will be used depending on the post-surgical phase, prosthetic phase or advanced phase.

Contraindications:

Placement of dental implants and prosthetic components is contraindicated in those patients who cannot be subjected to surgical interventions due to any medical pathology.

The MPI prosthetic components are single use devices, should not be re-used.
Sterile prosthetic components should not be re-sterilized. Infections and adverse reactions with the consequent non-oseointegration might be caused if a single use device, is re-used as the sterility levels offered by MPI are not guaranteed any more.

Contraindicated to patients who are being treated with radiotherapy, bisphosphonates and who have a diagnosed allergy to Titanium and the materials indicated in this IFU, infection or periodontal disease, limited bone availability, and smokers.

Patients, who had been subjected to radiotherapy or had been treated with bisphosphonates before the implantation procedure, need to contact the specialist who had treated them in order to check the risks and contraindications.

Related to patients diagnosed with bruxism, the professional should evaluate in his diagnosis and treatment plan the involuntary habit of clenching and grinding their teeth without any funtional purpose.
This should be reflected in the occlusal protection of the implant during the osseointegration phase and in the design of the definitive prosthesis.

Clinician must inform to the patients about contraindications and possible adverse reactions, as well as the importance of maintain a good dental hygiene.

Any serious incident occurring in connection with the device must be notified to the manufacturer and to the competent authority of the Member State in which the user and/or the patient is established.

Possible adverse reactions

- > Fenestrations and dehiscences
- > Infection
- > Loss of integration
- > Discomfort or pain
- > Inflammation
- > Bone loss.

To date there have been no known significant adverse effects caused by this medical device.

Warnings:

When placing a dental implant or a prosthetic component is necessary taking into account the medical situation of the patient, his physiological and anatomic conditions as well as his bone quality and quantity in the area.

It is highly recommended that the MPI implant system should only be placed by trained professionals.
The surgical and restorative techniques required to properly utilize these devices are highly specialized and complex procedures. Improper techniques can lead to non-osseointegration of the implant, bone loss and restoration fracture.

Even if the clinicians have obtained the necessary training in the surgical techniques of oral implantology, it is advisable that they consult the surgical and prosthetic protocols of MPI before using the MPI devices, as well as receiving personalized assistance by MPI.

Precautions:

During the diagnosis and planning phase, it is very important to take into account the indications and contraindications of this instruction and to know the MPI implant system and its prosthetic possibilities, in order to reduce as much as possible the potential error variables.

Once the clinical case has completely been planned and the implant's location, diameter and length been decided, the clinician will proceed by designing the type of prosthesis that will be used.

Subsequently, the patient will be operated, taking into account the incision design, the drilling protocol, the placement of the implant and the sterile conditions of the surgical area.

In case the device is not osseointegrated correctly, it can be safely disposed after use, clinical protocols should be followed to avoid infections and contaminations.

Sterilization, handling storage and transport:

The packaging has been validated simulating transport and storing conditions to guarantee protection against impacts and the sterility of the sterile devices.
Some prosthetic components are delivered in a sterile state (sterilized using gamma irradiation) and packed in a heat sealed blister.
Some prosthetic components are delivered non-sterile and packed in a heat sealed-pouches. The non-sterile devices are not intended to be used sterile.

These devices are packaged following concrete protocols to guarantee optimal results. Before using any product of the MPI implant system, the packaging should be examined looking out for possible damages that could affect the sterility of the product.

These devices should be used with the rest of articles and references of the MPI family in order to obtain optimal results.

MPI prosthetic components should be stored in original packaging, kept out of direct sunlight, avoid major thermal fluctuation and away from any source of artificial ultraviolet light.

Sterile devices packed in blister, should be stored in dry storage areas between 10-30°C and 50% humidity, throughout 5 years (expiration period to ensure its sterility).

Non-sterile devices packed in pouches, should be stored in dry storage areas between 15-25°C and 40-60% humidity, throughout 5 years (established period to guarantee the correct preservation of the device).

The product is perfectly packaged and heat sealed, a defect of this can mean the loss of decontamination and disinfection properties, so it is recommended to discard its use.

These are the ideal storing conditions advised for the MPI implant system.

For more information, please visit MPI's website: www.mpimplants.com.

Labelling symbols:

Manufacturer

Reference

Lot number

CE symbol with the notified body number 0120

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