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**1. Product:**  
MPI dental implant systems.

Implant systems types:

MPI Privilege® External Connection	
Reference	Description
IPHE310	MPI Privilege® HE dental implant Ø 3,3/3,5 x 10 mm
IPHE311	MPI Privilege® HE dental implant Ø 3,3/3,5 x 11,5 mm
IPHE313	MPI Privilege® HE dental implant Ø 3,3/3,5 x 13 mm
IPHE315	MPI Privilege® HE dental implant Ø 3,3/3,5 x 15 mm
IPHE3410	MPI Privilege® HE dental implant Ø 3,3/4,1 x 10 mm
IPHE3411	MPI Privilege® HE dental implant Ø 3,3/4,1 x 11,5 mm
IPHE3413	MPI Privilege® HE dental implant Ø 3,3/4,1 x 13 mm
IPHE3415	MPI Privilege® HE dental implant Ø 3,3/4,1 x 15 mm
IPHE47	MPI Privilege® HE dental implant Ø 4/4,1 x 7 mm
IPHE48	MPI Privilege® HE dental implant Ø 4/4,1 x 8,5 mm
IPHE410	MPI Privilege® HE dental implant Ø 4/4,1 x 10 mm
IPHE411	MPI Privilege® HE dental implant Ø 4/4,1 x 11,5 mm
IPHE413	MPI Privilege® HE dental implant Ø 4/4,1 x 13 mm
IPHE415	MPI Privilege® HE dental implant Ø 4/4,1 x 15 mm
IPHE57	MPI Privilege® HE dental implant Ø 5/5,1 x 7 mm

MPI Privilege® External Connection	
Reference	Description
IPHE510	MPI Privilege® HE dental implant Ø 5/5,1 x 10 mm
IPHE511	MPI Privilege® HE dental implant Ø 5/5,1 x 11,5 mm
IPHE513	MPI Privilege® HE dental implant Ø 5/5,1 x 13 mm
IPHE515	MPI Privilege® HE dental implant Ø 5/5,1 x 15 mm
IPHE58	MPI Privilege® HE dental implant Ø 5/5,1 x 8,5 mm

MPI Privilege® Internal Connection	
Reference	Description
IECM38	MPI Privilege® CM dental implant Ø 3,5 x 8,5 mm
IECM310	MPI Privilege® CM dental implant Ø 3,5 x 10 mm
IECM311	MPI Privilege® CM dental implant Ø 3,5 x 11,5 mm
IECM313	MPI Privilege® CM dental implant Ø 3,5 x 13 mm
IECM315	MPI Privilege® CM dental implant Ø 3,5 x 15 mm
IECM47	MPI Privilege® CM dental implant Ø 4 x 7 mm
IECM48	MPI Privilege® CM dental implant Ø 4 x 8,5 mm
IECM410	MPI Privilege® CM dental implant Ø 4 x 10 mm
IECM411	MPI Privilege® CM dental implant Ø 4 x 11,5 mm
IECM413	MPI Privilege® CM dental implant Ø 4 x 13 mm
IECM415	MPI Privilege® CM dental implant Ø 4 x 15 mm
IECM510	MPI Privilege® CM dental implant Ø 5 x 10 mm
IECM511	MPI Privilege® CM dental implant Ø 5 x 11,5 mm
IECM513	MPI Privilege® CM dental implant Ø 5 x 13 mm
IECM515	MPI Privilege® CM dental implant Ø 5 x 15 mm

MPI Excellence® Internal Connection	
Reference	Description
IECM38	MPI Excellence® CM dental implant Ø 3,5 x 8,5 mm
IECM310	MPI Excellence® CM dental implant Ø 3,5 x 10 mm
IECM311	MPI Excellence® CM dental implant Ø 3,5 x 11,5 mm
IECM313	MPI Excellence® CM dental implant Ø 3,5 x 13 mm
IECM315	MPI Excellence® CM dental implant Ø 3,5 x 15 mm
IECM378	MPI Excellence® CM dental implant Ø 3,75 x 8,5 mm
IECM3710	MPI Excellence® CM dental implant Ø 3,75 x 10 mm
IECM3711	MPI Excellence® CM dental implant Ø 3,75 x 11,5 mm
IECM3713	MPI Excellence® CM dental implant Ø 3,75 x 13 mm
IECM3715	MPI Excellence® CM dental implant Ø 3,75 x 15 mm
IECM48	MPI Excellence® CM dental implant Ø 4,3 x 8,5 mm
IECM410	MPI Excellence® CM dental implant Ø 4,3 x 10 mm
IECM411	MPI Excellence® CM dental implant Ø 4,3 x 11,5 mm
IECM413	MPI Excellence® CM dental implant Ø 4,3 x 13 mm
IECM415	MPI Excellence® CM dental implant Ø 4,3 x 15 mm
IECM510	MPI Excellence® CM dental implant Ø 5 x 10 mm
IECM511	MPI Excellence® CM dental implant Ø 5 x 11,5 mm
IECM513	MPI Excellence® CM dental implant Ø 5 x 13 mm
IECM515	MPI Excellence® CM dental implant Ø 5 x 15 mm

MPI Short® Internal Connection	
Reference	Description
ISCM56	MPI Short® CM dental implant Ø 5 x 6 mm
ISCM57	MPI Short® CM dental implant Ø 5 x 7 mm
ISCM58	MPI Short® CM dental implant Ø 5 x 8,5 mm

MPI All-In® Internal Connection	
Reference	Description
IACM38	MPI All-In® CM dental implant Ø 3,5 x 8,5 mm
IACM310	MPI All-In® CM dental implant Ø 3,5 x 10 mm
IACM311	MPI All-In® CM dental implant Ø 3,5 x 11,5 mm
IACM313	MPI All-In® CM dental implant Ø 3,5 x 13 mm
IACM315	MPI All-In® CM dental implant Ø 3,5 x 15 mm
IACM378	MPI All-In® CM dental implant Ø 3,75 x 8,5 mm
IACM3710	MPI All-In® CM dental implant Ø 3,75 x 10 mm
IACM3711	MPI All-In® CM dental implant Ø 3,75 x 11,5 mm
IACM3713	MPI All-In® CM dental implant Ø 3,75 x 13 mm
IACM3715	MPI All-In® CM dental implant Ø 3,75 x 15 mm
IACM48	MPI All-In® CM dental implant Ø 4,3 x 8,5 mm
IACM410	MPI All-In® CM dental implant Ø 4,3 x 10 mm
IACM411	MPI All-In® CM dental implant Ø 4,3 x 11,5 mm
IACM413	MPI All-In® CM dental implant Ø 4,3 x 13 mm
IACM415	MPI All-In® CM dental implant Ø 4,3 x 15 mm
IACM510	MPI All-In® CM dental implant Ø 5 x 10 mm
IACM511	MPI All-In® CM dental implant Ø 5 x 11,5 mm
IACM513	MPI All-In® CM dental implant Ø 5 x 13 mm

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Dental implants are medical devices (dental), subgingival endosseous implants made of commercially pure Titanium Grade IV cold drawn bars.

Dental implants are medical devices class IIb.

## 3. Intended use / Purpose:

Dental implants are designed and conceived to be surgically implanted in the maxilla or jaw, in order to provide fixation to support and retain dental prostheses to achieve an aesthetic and functional improvement, a better distribution of chewing forces and a decrease in bone loss.

</

necessary protective and clothing measures. Improper use can lead to injure the patient.

- They are intended to be used in a sterile environment.
- The dental implants indicated must be used with the rest of the articles and references of the system to which they belong in order to guarantee a good result.
- Intended users should review the expiration date indicated on the package and check that the package is not damaged or opened to maintain the sterility.
- Dental implants are available in different lengths and diameters in order to be able to rehabilitate different cases according to vertical and horizontal bone availability.
- The implants are delivered packaged with their corresponding cover screw (included in the tube cap). Before placing the implant, the intended user should verify that the cover screw is suitable for that implant, and the corresponding constants for them.
- MPI Short® dental implants are recommended for multiple restorations, they are not recommended to immediate loading.

- When it is decided to discard an implant, the intended user must dispose it safely following the protocol for the elimination of medical devices defined in the dental clinic, in accordance with hospital regulations and local existing legislation, to guarantee the safety elimination and avoid infections and contamination.

- Maximum recommended torques depending on the dental implant type (established according to the mechanical tests carried out):

Type of implant	Maximum recommended torques (N·cm)
MPI Privilege® Internal Connection	80
MPI Privilege® External Connection	110
MPI Excellence® Internal Connection	80
MPI Short® Internal Connection	80
MPI All-in® Internal Connection	80

[\*]

#### 8. Contraindications:

It is contraindicated to use the dental implants in patients:
- Not suitable for an oral surgical procedure (for example, patients who have undergone radiotherapy or who have been treated with bisphosphonates, should contact the specialist who has treated them to verify the risks and contraindications).

- Allergic or hypersensitive to the materials used.

- With periodontal diseases, limited bone availability and density, diabetics and/or smokers and/or bruxism. In case of a patient with bruxism (the involuntary habit of clenching and grinding the teeth without functional purpose), the intended user should assess it in the diagnosis and also in the treatment plan because it should be reflected in the occlusal protection of the implant during bone integration, and in the design of the final prosthesis.

- The placement of dental implants is contraindicated in those patients who cannot undergo surgery for any medical condition.

- The intended user must inform the intended patient of these contraindications.

##### 9. Possible adverse reactions and risks:

- An improper use of the products, an inadequate drilling protocol/surgical technique, can cause bone loss, fenestrations, dehiscences, bone and/or nerve damage, discomfort and pain in the patient, and the consequent non-osseointegration of the implant.

- An inadequate drilling protocol with insufficient irrigation (not cooling the instruments with sterile saline solution) could also cause necrosis in the bone due to thermal effects.

- The inadequate application of higher torque values than the recommended ones, could damage the connection between the implant and the contra-angle; in addition, it could lead to a lack of primary stability, and the bone could suffer necrosis, which would mean failure in the implant osseointegration.

- The unskilled or untrained user, without planning clinical cases and/or in contraindicated patients, may cause harm to the patient.

- An inappropriate environment where surgery is to be performed/improper storage conditions, can cause contamination of the product, and can lead to infection of the patient and/or non-osseointegration of the implant.

To date, no other possible significant adverse effects caused by these products are known.

The intended user must inform the intended patient about the importance of maintaining good dental hygiene, and about these possible adverse reactions and risks.

##### 10. Sterilization:

The indicated dental implants are delivered sterile. They are sterilized by gamma irradiation and they are intended for single use.

##### 11. Precautions:

- Before placing the dental implant, check that it is sterilized (the irradiation label must be read).

- It is recommended that the intended user has received regulated and supervised training by expert professionals on surgical techniques in oral implantology; the intended user should consult the specific surgical and prosthetic protocols; the intended user should receive personalized assistance and the intended user should attend special training before using the products. MPI offers courses for different levels of knowledge and experience.

- The intended user who uses these products must have experience to use them safely and properly according to these instructions for use, in addition the intended user must inform intended patients about contraindications and possible adverse reactions, as well as the importance of maintaining good dental hygiene.

- During the phase of the diagnosis and planning the implant case, the indications and contraindications of these instructions for use will be considered, the MPI implant system and its prosthetic possibilities should be known, in order to minimize error variables.

- Once the complete case has been planned, and the location, the diameter and the length of the implant have been decided, the type of prosthesis to be used will be designed.

- Subsequently, the patient will undergo surgery, taking into account the design of the incision, the drilling protocol, the placement of the implant and the sterilization conditions of the operative field.

[\*]

##### 12. Warnings:

- It is necessary to consider the medical situation of the patient, the physiological and anatomical conditions, and the bone quality and bone quantity in the area where the implant is going to be implanted. An inadequate technique can incur any of the possible adverse reactions and risks mentioned above.

- This product is for single use, it cannot be restertilized or re-used. If dental implants are re-used, the sterilization levels offered by MPI are not guaranteed, as they can cause infections and adverse reactions in the patient and therefore the non-osseointegration of the implant.

- The product is perfectly packaged and heat-sealed. Check if the sealed package is damaged or opened because the sterility and cleaning conditions of the product could have been altered; In that case, do not use it.

- All products have a batch number to allow their traceability, providing greater security to the intended user, due to with this batch number it is possible to know the entire history of the product from the manufacturing process to the moment of its distribution.

- This product is intended to be used in a sterile environment.

- Any serious incident related to the device must be notified to the manufacturer and to the competent authority of the Member State in which the intended user and/or the patient is established.

##### 13. Handling:

- To place the implants, the implant motor connector must be used; to activate the rotation of this instrument, they must be connected to the handpiece (contra-angle) of the implant motor.

- To accurately measure the drilling mode with the dental drills, it is recommended to verify the drilling depth using calibrated depth gauges.

- To avoid mechanical and thermal risks, use abundant irrigation with sterile saline solution continuously throughout all the implant insertion process.

##### 14. Storage and transport:

The implants are wrapped and packaged according to defined and validated protocols that guarantee optimal results. The original packaging has been validated simulating the transport and storage conditions to guarantee the cleanliness of the packaged products, their protection against impacts, and the sterility of the implants. Recommended ideal storage conditions: MPI dental implants should be stored in its original package in a clean and dry storage areas between 10°C-30°C and 50% humidity, avoiding temperature fluctuations, not exposed to direct sunlight and keep away from any source of artificial ultraviolet light for 5 years (expiration period to ensure its sterility).

Incorrect storage during transport can influence the characteristics of the device and lead to its failure.

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